ANTIBIOGRAMMA 2012

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1 CLSI Performance Standards for Antimicrobial Susceptibility Testing. 20th Informational Supplement - M100 S20. January 2010.

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"-" indicates that susceptibility testing is not recommended as the species is a poor target for therapy with the drug. Isolates may be reported as R without prior testing.

"IE" indicates that there is insufficient evidence that the species in question is a good target for therapy with the drug. An MIC with a comment but without an accompanying S, I or R categorisation may be reported

ENTEROBACTERIACEAE	EUCAST		CLSI	
PENICILLINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ampicillin	8	8	8	32
Amoxicillin-clavulanate	8	8	8/4	32/16
Piperacillin-tazobactam	8	16	16/4	128/4

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefepime	1	4	8	32
Cefotaxime	1	2	1	4
Ceftazidime	1	4	4	16
Ceftriaxone	1	2	1	4

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	0.5	1	0.5	2
Imipenem	2	8	1	4
Meropenem	2	8	1	4

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	0.5	1	1	4
Levofloxacin	1	2	2	8
Norfloxacin	0.5	1	4	16


AMINOGLYCOSIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Amikacin	8	16	16	64
Gentamicin	2	4	4	16
Tobramycin	2	4	4	16

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	-	-	4	16
Tigecycline ¹	1	2	-	-

MISCELLANEOUS AGENTS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Colistin	2	2	-	-
Fosfomycin iv	32	32	64	256
Nitrofurantoin (uncomplicated UTI only)	64 ²	64 ²	32	128
Trimethoprim-sulfamethoxazole ²	2	4	2/38	4/76

Note 1 - Tigecycline has decreased activity against *Morganella* spp and *Providencia* spp.

Note 2 - Breakpoint relate to *E.coli* only

PSEUDOMONAS SPP.	EUCAST	CLSI
-------------------------	---------------	-------------

PENICILLINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Piperacillin	16	16	16	128
Piperacillin-tazobactam	16	16	16/4	128/4

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefepime	8	8	8	32
Ceftazidime	8	8	8	32

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	-	-	-	-
Imipenem	4	8	2	8
Meropenem	2	8	2	8

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	0.5	1	1	4
Levofloxacin	1	2	2	8

AMINOGLYCOSIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Amikacin	8	16	16	64
Gentamicin	4	4	4	16
Tobramycin	4	4	4	16

MISCELLANEOUS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Colistin	4	4	2	8

STENOTROPHOMONAS MALTOPHILIA	EUCAST	CLSI
-----------------------------------------	---------------	-------------

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Levofloxacin	-	-	2	8

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Minocycline	-	-	4	16

MISCELLANEOUS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Trimethoprim-sulfamethoxazole	4	4	2/38	4/76

ACINETOBACTER SPP.	EUCAST	CLSI
--------------------	--------	------

PENICILLINS ¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ampicillin-sulbactam	IE	IE	8/4	32/16
Piperacillin	IE	IE	16	128
Piperacillin-tazobactam	IE	IE	16/4	128/4

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	-	-	-	-
Imipenem	2	8	4	16
Meropenem	2	8	4	16

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	1	1	1	4
Levofloxacin	1	2	2	8

AMINOGLYCOSIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Amikacin	8	16	16	64
Gentamicin	4	4	4	16
Tobramycin	4	4	4	16

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tigecycline	IE	IE	-	-

MISCELLANEOUS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Colistin	2	2	2	4
Rifampicin	-	-	-	-
Trimethoprim-sulfamethoxazole	2	4	2/38	4/76

Note 1: Susceptibility testing of *Acinetobacter* spp. to penicillins is unreliable. In most instances *Acinetobacter* spp. are resistant to penicillins.

STAPHYLOCOCCUS SPP.	EUCAST	CLSI
--------------------------------	---------------	-------------

PENICILLINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin	0.12 ¹	0.12 ¹	0.12 ^A	0.25 ^A
Ampicillin	Note ¹	Note ¹	0.25	0.5
Ampicillin-sulbactam	Note ¹	Note ¹	8/4	32/16
Amoxicillin	Note ¹	Note ¹	Note ^A	Note ^A
Amoxicillin-clavulanate	Note ¹	Note ¹	4/2	8/4
Piperacillin	Note ¹	Note ¹	Note ^A	Note ^A
Piperacillin-tazobactam	Note ¹	Note ¹	8/4	16/4
Oxacillin²	Note ²	Note ²	0.25	0.5

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefepime	Note ³	Note ³	8	32
Cefotaxime	Note ³	Note ³	8	64
Cefoxitin (screen)³	Note ³	Note ³	8	16

CARBAPENEMS ¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	Note ¹	Note ¹	2	8
Imipenem	Note ¹	Note ¹	4	16
Meropenem	Note ¹	Note ¹	4	16

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	1	1	1	4
Levofloxacin	1	2	1	4

AMINOGLYCOSIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Amikacin	8	16	16	64
Gentamicin	1	1	4	16
Tobramycin	1	1	4	16

GLYCOPEPTIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Teicoplanin, <i>S. aureus</i>	2	2	8	32
Teicoplanin, Coagulase-negative staphylococci	4	4	8	32
Vancomycin, <i>S. aureus</i>	2	2	2	16
Vancomycin, Coagulase-negative staphylococci	4	4	4	32

MACROLIDES AND LINCOSAMIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Erythromycin ⁴	1	2	0.5	8
Clindamycin	0.25	0.5	0.5	4

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	1	2	4	16
Tigecycline	0.5	0.5	-	-

MISCELLANEOUS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Daptomycin	1	1	1	
Fusidic acid	1	1	-	-
Linezolid	4	4	4	8
Rifampicin	0.06	0.5	1	4
Trimethoprim-sulfamethoxazole	2	4	2/38	4/76

Note 1: The beta-lactam (with or without inhibitor) susceptibility is inferred from the oxacillin susceptibility.

Note 2: *S. aureus* and *S. lugdunensis* with oxacillin MIC values >2 mg/L are mostly methicillin resistant due to the presence of the *mecA* gene. The corresponding oxacillin MIC for coagulase-negative staphylococci is >0.25 mg/L.

Note 3: Susceptibility of staphylococci to cephalosporins is inferred from the cefoxitin susceptibility except for ceftazidime, which do not have breakpoints and should not be used for staphylococcal infections. *S. aureus* and *S. lugdunensis* with cefoxitin MIC values >4 mg/L are mostly methicillin resistant due to the presence of the *mecA* gene. For coagulase-negative staphylococci other than *S. lugdunensis* the cefoxitin MIC is a poorer predictor of methicillin resistance than the disk diffusion test.

Note 4: Erythromycin can be used to determine susceptibility to azithromycin, clarithromycin and roxithromycin.

Note A: Penicillin should be used to test the susceptibility of all penicillinase-labile penicillins such as amoxicillin, piperacillin

ENTEROCOCCUS SPP.	EUCAST	CLSI
------------------------------	---------------	-------------

IN ENDOCARDITIS, REFER TO NATIONAL OR INTERNATIONAL
ENDOCARDITIS GUIDELINES FOR BREAKPOINTS FOR ENTEROCOCCUS
SPP

PENICILLINS¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin	-	-	-	-
Ampicillin	4	8	8	16
Ampicillin-sulbactam²	4	8	Note ^A	Note ^A

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefepime	-	-	-	-
Cefotaxime	-	-	-	-
Ceftazidime	-	-	-	-
Ceftriaxone	-	-	-	-

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	-	-	-	-
Imipenem	4	8	Note ^A	Note ^A
Meropenem	-	-	-	-

AMINOGLYCOSIDES ³	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Gentamicin	Note ⁴	Note ⁴	-	-
Tobramycin	IE	IE	-	-

GLYCOPEPTIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Teicoplanin	2	2	8	32
Vancomycin	4	4	4	32

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	-	-	4	16
Tigecycline	0.25	0.5	-	-

MISCELLANEOUS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Linezolid	4	4	2	8
Nitrofurantoin (uncomplicated UTI only)	64 ⁵	64 ⁵	32	128
Trimethoprim-sulfamethoxazole	0.03	1	-	-

Note 1: *E. faecium* resistant to penicillins can be considered resistant to all other beta-lactam agents including carbapenems.

Note 2: Susceptibility to ampicillin, amoxicillin and piperacillin with and without beta-lactamase inhibitor can be inferred from the ampicillin susceptibility test.

Note 3: Aminoglycoside monotherapy is ineffective against enterococci. There is synergism between aminoglycosides and beta-lactam agents against enterococci without acquired aminoglycoside resistance mechanisms.

Note 4: Isolates with gentamicin MIC >128 mg/L have acquired resistance mechanisms and can be reported as high-level aminoglycoside resistant (with the exception of streptomycin, which must be tested separately). There is no synergistic effect between aminoglycosides and beta-lactam agents in enterococci with high-level aminoglycoside resistance.

Note 5: Nitrofurantoin breakpoints are valid for *E. faecalis* only.

Note A: Ampicillin is the class representative for ampicillin and amoxicillin. Ampicillin result may be used to predict susceptibility to amoxicillin-clavulanic acid, ampicillin-sulbactam, piperacillin and piperacillin-tazobactam among non β -lactamase producing enterococci. Ampicillin can be used to predict imipenem susceptibility providing the species is confirmed to be *E. faecalis*.

STREPTOCOCCUS GROUPS A, B, C AND G	EUCAST	CLSI
---------------------------------------------------	---------------	-------------

PENICILLINS¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin	0.25	0.25	0.12	
Ampicillin	Note ¹	Note ¹	0.25	0.25
Ampicillin-sulbactam²	Note ¹	Note ¹	Note ^A	Note ^A
Amoxicillin	Note ¹	Note ¹	Note ^A	Note ^A
Amoxicillin-clavulanate²	Note ¹	Note ¹	Note ^A	Note ^A

CEPHALOSPORINS¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefotaxime	Note ¹	Note ¹	Note ^A	Note ^A
Ceftazidime	-	-	-	-
Ceftriaxone	Note ¹	Note ¹	Note ^A	Note ^A
Cefuroxime	Note ¹	Note ¹	Note ^A	Note ^A

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	-	-	-	-
Levofloxacin	1	1	2	8

GLYCOPEPTIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Teicoplanin	2	2	2	
Vancomycin	2	2	1	

MACROLIDES AND LINCOSAMIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Erythromycin ³	0.25	0.25	0.25	1
Clindamycin	0.5	0.5	0.25	1

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	1	1	2	8

MISCELLANEOUS AGENTS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Linezolid	2	2	2	
Nitrofurantoin (uncomplicated UTI only)	64 ⁴	64 ⁴	-	-

Note 1: The beta-lactam susceptibility of beta-haemolytic streptococcus groups A, B, C and G is inferred from the penicillin susceptibility.

Note 2: Streptococcus groups A, B, C and G do not produce beta-lactamase. The addition of a beta-lactamase inhibitor does not add clinical benefit.

Note 3: Erythromycin can be used to determine susceptibility to azithromycin, clarithromycin and roxithromycin.

Note 4: Nitrofurantoin breakpoints apply to *S. agalactiae* (group B streptococci) only.

Note A: An organism susceptible to penicillin can be considered susceptible to Ampicillin amoxicillin-clavulanic acid, ampicillin-sulbactam, piperacillin and piperacillin-tazobactam, cefotaxime, ceftazidime, ceftriaxone e cefuroxime.

STREPTOCOCCUS PNEUMONIAE	EUCAST	CLSI
-------------------------------------	---------------	-------------

PENICILLINS¹	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin (infections other than meningitis)	0.06 ²	2 ²	2	8
Benzylpenicillin (meningitis)	0.06	0.06	0.06	0.12
Ampicillin	0.5 ¹	2 ¹	-	-
Amoxicillin	Note ¹	Note ¹	2 ^A	8 ^A
Amoxicillin-clavulanate	Note ¹	Note ¹	2/1 ^A	8/4 ^A
Piperacillin	Note ¹	Note ¹	-	-
Piperacillin-tazobactam	Note ¹	Note ¹	-	-

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefotaxime	0.5 ¹	2	1 ^A /0.5 ^B	4 ^A /2 ^B
Ceftazidime	-	-	-	-
Ceftriaxone	0.5 ¹	2	1 ^A /0.5 ^B	4 ^A /2 ^B

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem ³	0.5	0.5	1	4
Imipenem ³	2	2	0.12	1
Meropenem (infections other than meningitis)	2	2	0.25	1
Meropenem (meningitis)	0.25	1	0.25	1

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	0.12	2	-	-
Levofloxacin ⁴	2	2	2	8

GLYCOPEPTIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Teicoplanin	2	2	-	-
Vancomycin	2	2	1	

MACROLIDES AND LINCOSAMIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Erythromycin ⁵	0.25	0.5	0.25	1
Clindamycin	0.5	0.5	0.25	1

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	1	2	2	8

MISCELLANEOUS AGENTS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Trimethoprim-sulfamethoxazole	1	2	0.5/9.5	4/76

Note 1: Most MIC values for penicillin, ampicillin, amoxicillin and piperacillin (with or without a beta-lactamase inhibitor) do not significantly differ and can be reported susceptible to beta-lactam agents that have been given breakpoints.

Note 2: In pneumonia, when a dose of 1.2 g x 4 is used, isolates with MIC ≤0.5 mg/L should be regarded as susceptible to benzylpenicillin. In pneumonia, when a dose of 2.4 g x 4 or 1.2 g x 6 is used, isolates with MIC ≤1 mg/L should be regarded as susceptible to benzylpenicillin. In pneumonia, when a dose of 2.4 g x 6 is used, isolates with MIC ≤2 mg/L should be regarded as susceptible.

Note 3. Not for meningitis (meropenem is the only carbapenem used for meningitis).

Note 4: The breakpoints for levofloxacin relate to high dose therapy.

Note 5: Erythromycin can be used to determine susceptibility to azithromycin, clarithromycin and roxithromycin.

Note A: non- meningitis

Note B: meningitis

VIRIDANS GROUP STREPTOCOCCI	EUCAST	CLSI
----------------------------------------	---------------	-------------

IN ENDOCARDITIS, REFER TO NATIONAL OR INTERNATIONAL ENDOCARDITIS GUIDELINES FOR BREAKPOINTS FOR VIRIDANS GROUP STREPTOCOCCI.

PENICILLINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin	0.25	2	0.12	4
Ampicillin	0.5	2	0.25	8
Amoxicillin-clavulanate	Note ¹	Note ¹	-	-

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefotaxime	0.5	0.5	1	4
Ceftriaxone	0.5	0.5	1	4

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem	0.5	0.5	1	
Imipenem	2	2	-	-
Meropenem	2	2	0.5	

GLYCOPEPTIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Teicoplanin	2	2	-	-
Vancomycin	2	2	1	

MACROLIDES AND LINCOSAMIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Erythromycin	IE	IE	0.25	1
Clindamycin	0.5	0.5	0.25	1

Note 1: For isolates susceptible to benzylpenicillin, susceptibility can be inferred from benzylpenicillin or ampicillin. For isolates resistant to benzylpenicillin, susceptibility is inferred from ampicillin.

HAEMOPHILUS INFLUENZAE	EUCAST	CLSI
-----------------------------------	---------------	-------------

PENICILLINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Benzylpenicillin	IE	IE	IE	IE
Ampicillin	1	1	1	4
Amoxicillin	2	2	Note ^A	Note ^A
Amoxicillin-clavulanate	2	2	4/2	8/4

CEPHALOSPORINS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Cefotaxime	0.12	0.12	2	
Ceftazidime	-	-	2	
Ceftriaxone	0.12	0.12	2	

CARBAPENEMS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ertapenem¹	0.5	0.5	0.5	
Imipenem¹	2	2	4	
Meropenem (infections other than meningitis)	2	2	0.5	
Meropenem (meningitis)	0.25	1	0.5	

FLUOROQUINOLONES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Ciprofloxacin	0.5	0.5	1	
Levofloxacin	1	1	2	

MACROLIDES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Erythromycin	0.5	16		

TETRACYCLINES	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Tetracycline	1	2	2	8

MISCELLANEOUS AGENTS	MIC breakpoint (mg/L)		MIC breakpoint (mg/L)	
	S ≤	R >	S ≤	R ≥
Rifampicin (for prophylaxis only)	1	1	1	4
Trimethoprim-sulfamethoxazole	0.5	1	0.5/9.5	4/76

Note 1: Not for meningitis (meropenem is the only carbapenem used for meningitis).

Note A: The results of ampicillin susceptibility tests should be used to predict the activity of amoxicillin.

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REVIEW
EUCAST expert rules antimicrobial susceptibility testing

R.Leclercq^{1,2}, R. Cantòn^{2,3,4}, D.F.J. Brown⁴, C.G. Giske^{2,4,5}, P. Heising^{2,4}, A.P. Mac Gowan^{4,7}, J. W. Mouton^{4,8}, P. Nordmann^{2,9}, A. C. Rodloff^{4,10}, G.M. Rossolini^{2,11}, C.J. Soussy^{4,12}, M. Steinbakk^{4,13}, T.G. Winstanley^{2,14}, and G. Kahlmeter^{4,15}

Clin Microbiol Infect. 2011 Oct 21. doi: 10.1111/j.1469-0691.2011.03703.x.
<http://onlinelibrary.wiley.com/doi/10.1111/j.1469-0691.2011.03703.x/full>

TABLE 1. Intrinsic resistance in Enterobacteriaceae; Enterobacteriaceae are also intrinsically resistant to benzylpenicillin, glycopeptides, fusidic acid, macrolides (with some exceptions a), lincosamides, streptogramins, rifampicin, daptomycin, and linezolid

Rule no.	Organisms	Ampicillin	Amoxycillin clavulanate	Ticarcillin
1.01	Citrobacter koseri	R	–	R
1.02	Citrobacter freundii	R	R	–
1.03	Enterobacter cloacae	R	R	–
1.04	Enterobacter aerogenes	R	R	–
1.05	Escherichia hermannii	R	–	R
1.06	Hafnia alvei	R	R	–
1.07	Klebsiella spp.	R	–	R
1.08	Morganella morganii	R	R	–
1.09	Proteus mirabilis	–	–	–
1.10	Proteus vulgaris	R	–	–
1.11	Proteus penneri	R	–	–
1.12	Providencia rettgeri	R	R	–
1.13	Providencia stuartii	R	R	–
1.14	Serratia marcescens	R	R	–
1.15	Yersinia enterocolitica	R	R	R
1.16	Yersinia pseudotuberculosis	–	–	–

Rule no.	Organisms	Piperacillin	Cefazolin	Cefoxitin
1.01	Citrobacter koseri	R	–	–
1.02	Citrobacter freundii	–	R	R
1.03	Enterobacter cloacae	–	R	R
1.04	Enterobacter aerogenes	–	R	R
1.05	Escherichia hermannii	–	–	–
1.06	Hafnia alvei	–	R	–
1.07	Klebsiella spp.	–	–	–
1.08	Morganella morganii	–	R	–
1.09	Proteus mirabilis	–	–	–
1.10	Proteus vulgaris	–	R	–
1.11	Proteus penneri	–	R	–
1.12	Providencia rettgeri	–	R	–
1.13	Providencia stuartii	–	R	–
1.14	Serratia marcescens	–	R	–
1.15	Yersinia enterocolitica	–	R	R
1.16	Yersinia pseudotuberculosis	–	–	–

Rule no.	Organisms	Cefamandole	Cefuroxime	Aminoglycosides
1.01	Citrobacter koseri	–	–	–
1.02	Citrobacter freundii	–	–	–
1.03	Enterobacter cloacae	–	–	–
1.04	Enterobacter aerogenes	–	–	–
1.05	Escherichia hermannii	–	–	–
1.06	Hafnia alvei	–	–	–
1.07	Klebsiella spp.	–	–	–
1.08	Morganella morganii	–	R	–
1.09	Proteus mirabilis	–	–	–
1.10	Proteus vulgaris	R	R	–
1.11	Proteus penneri	R	R	–
1.12	Providencia rettgeri	–	–	–
1.13	Providencia stuartii	–	–	Noteb
1.14	Serratia marcescens	R	R	Notec
1.15	Yersinia enterocolitica	R	–	–
1.16	Yersinia pseudotuberculosis	–	–	–

Rule no.	Organisms	Tetracyclines/ tigecycline	PolymyxinB colistin	Nitrofurantoin
1.01	Citrobacter koseri	–	–	–
1.02	Citrobacter freundii	–	–	–
1.03	Enterobacter cloacae	–	–	–
1.04	Enterobacter aerogenes	–	–	–
1.05	Escherichia hermannii	–	–	–
1.06	Hafnia alvei	–	–	–
1.07	Klebsiella spp.	–	–	–
1.08	Morganella morganii	R	R	R
1.09	Proteus mirabilis	R	R	R
1.10	Proteus vulgaris	R	R	R
1.11	Proteus penneri	R	R	R
1.12	Providencia rettgeri	R	R	R
1.13	Providencia stuartii	R	R	R
1.14	Serratia marcescens	–	R	R
1.15	Yersinia enterocolitica	–	–	–
1.16	Yersinia pseudotuberculosis	–	R	–

R. resistant.

a. Azithromycin is effective in vivo for the treatment of typhoid fever, and erythromycin may be used to treat travellers' diarrhoea.

b. *Providencia stuartii* produces a chromosomal AAC(2')-Ia enzyme and should be considered to be resistant to clinically available aminoglycosides, except amikacin, arbekacin, and streptomycin. Some isolates express the enzyme poorly and can appear to be susceptible to netilmicin in vitro, but should be reported as resistant, as mutation can result in overproduction of this enzyme.

c. All *Serratia marcescens* isolates produce a chromosomal AAC(6')-Ic enzyme that affects the activity of clinically available aminoglycosides, except streptomycin, gentamicin, and arbekacin.

TABLE 2. Intrinsic resistance in non-fermentative Gram-negative bacteria; non-fermentative Gram-negative bacteria are also intrinsically resistant to benzylpenicillin, cefoxitin, cefamandole, cefuroxime, glycopeptides, fusidic acid, macrolides, lincosamides, streptogramins, rifampicin, daptomycin, and linezolid

Rule no.	Organisms	Ampicillin	Amoxicillin clavulanate	Ticarcillin	Ticarcillin clavulanate	Piperacillin	Piperacillin tazobactam	Cefazolin
2.01	Acinetobacter baumannii,	Ra	Ra	–	–	–	–	R
	Acinetobacter calcoaceticus							
2.02	Achromobacter xylosoxidans	R	–	–	–	–	–	R
2.03	Burkholderia cepacia complex^b	R	R	R	R	–	–	R
2.04	Elizabethkingia meningoseptica	R	–	R	R	–	–	R
2.05	Ochrobactrum anthropi	R	R	R	R	R	R	R
2.06	Pseudomonas aeruginosa	R	R	–	–	–	–	R
2.07	Stenotrophomonas maltophilia	R	R	R	–	R	R	R

Rule no.	Organisms	Cefotaxime	Ceftriaxone	Ceftazidime	Ertapenem	Imipenem	Meropenem	Ciprofloxacin
2.01	Acinetobacter baumannii,	R	R	-	R	-	-	-
	Acinetobacter calcoaceticus							
2.02	Achromobacter xylosoxidans	R	R	-	R	-	-	-
2.03	Burkholderia cepacia complex^b	-	-	-	R	R	-	R
2.04	Elizabethkingia meningoseptica	R	R	R	R	R	R	-
2.05	Ochrobactrum anthropi	R	R	R	R	-	-	-
2.06	Pseudomonas aeruginosa	R	R	-	R	-	-	-
2.07	Stenotrophomonas maltophilia	R	R	Rf	R	R	R	-

Rule no.	Organisms	Chloramphenicol	Aminoglycosides	Trimethoprim	Trimethoprim-sulfamethoxazole	Fosfomycin	Tetracyclines tigecycline	Polymyxin	B/colistin
2.01	Acinetobacter baumannii,	-	-	R	-	R	-	-	
	Acinetobacter calcoaceticus								
2.02	Achromobacter xylosoxidans	-	-	-	-	-	-	-	
2.03	Burkholderia cepacia complexb	R	Rc	R	-	R	-	R	
2.04	Elizabethkingia meningoseptica	-	-	-	-	-	-	R	
2.05	Ochrobactrum anthropi	-	-	-	-	-	-	-	
2.06	Pseudomonas aeruginosa	R	Noted	Re	Re	-	R	-	
2.07	Stenotrophomonas maltophilia	-	Rc	Rg	-	R	-	-	

R. resistant.

a. *Acinetobacter baumannii* may appear to be susceptible to ampicillin–sulbactam, owing to the activity of sulbactam against this species.

b. *Burkholderia cepacia* complex includes different species. Some strains may appear to be susceptible to some b-lactams in vitro, but they are clinically resistant and are shown as R in the table.

c. *Burkholderia cepacia* and *Stenotrophomonas maltophilia* are intrinsically resistant to all aminoglycosides. Intrinsic resistance is attributed to poor permeability and putative efflux. In addition, most *Stenotrophomonas maltophilia* isolates produce the AAC(6′)-Iz enzyme.

d. *Pseudomonas aeruginosa* is intrinsically resistant to kanamycin and neomycin, owing to low-level APH(3′)-IIb activity.

e. *Pseudomonas aeruginosa* is typically resistant to trimethoprim and moderately susceptible to sulfonamides. Although it may appear to be susceptible in vitro to trimethoprim–sulphamethoxazole, it should be considered to be resistant.

f. *Stenotrophomonas maltophilia* may show low ceftazidime MIC values but should be considered to be resistant.

g. *Stenotrophomonas maltophilia* is typically susceptible to trimethoprim–sulphamethoxazole but resistant to trimethoprim alone.

TABLE 3. Intrinsic resistance in Gram-negative bacteria other than Enterobacteriaceae and non-fermentative Gram-negative bacteria; Gram-negative bacteria other than Enterobacteriaceae and non-fermentative Gram-negative bacteria listed are also intrinsically resistant to glycopeptides, lincosamides, daptomycin, and linezolid

Rule no.	Organisms	Macrolides	Fusidic acid	Streptogramins	Trimethoprim	Nalidixic acid
3.01	Haemophilus influenzae	I	R	–	–	–
3.02	Moraxella catarrhalis	–	–	–	R	–
3.03	Neisseria spp.	–	–	–	R	–
3.04	Campylobacter fetus	–	R	R	R	R
3.05	Campylobacter jejuni Campylobacter coli	–	R	R	R	–

R. resistant; **I.** intermediate.

TABLE 4. Intrinsic resistance in Gram-positive bacteria; Gram-positive bacteria are also intrinsically resistant to aztreonam, temocillin, polymyxin B/colistin, and nalidixic Acid

Rule no.	Organisms	Fusidic acid	Ceftazidime	Cephalosporins
				(except ceftazidime)
4.01	Staphylococcus saprophyticus	R	R	–
4.02	Staphylococcus cohnii,	–	R	–
	Staphylococcus xylosus			
4.03	Staphylococcus capitis	–	R	–
4.04	Other coagulase-negative staphylococci and Staphylococcus aureus	–	R	–
4.05	Streptococcus spp.	R	–	–
4.06	Enterococcus faecalis	R	R	R
4.07	Enterococcus gallinarum, Enterococcus casseliflavus	R	R	R
4.08	Enterococcus faecium	R	R	R
4.09	Corynebacterium spp.	–	–	–
4.10	Listeria monocytogenes	–	R	R
4.11	Leuconostoc spp., Pediococcus spp.	–	–	–
4.12	Lactobacillus spp. (some species)	–	–	–
4.13	Clostridium ramosum, Clostridium innocuum	–	–	–

Rule no.	Organisms	Aminoglycosides	Erythromycin	Clindamycin
4.01	Staphylococcus saprophyticus	–	–	–
4.02	Staphylococcus cohnii,	–	–	–
	Staphylococcus xylosus			
4.03	Staphylococcus capitis	–	–	–
4.04	Other coagulase-negative staphylococci and Staphylococcus aureus	–	–	–
4.05	Streptococcus spp.	Ra	–	–
4.06	Enterococcus faecalis	Ra	R	R
4.07	Enterococcus gallinarum, Enterococcus casseliflavus	Ra	R	R
4.08	Enterococcus faecium	Ra,b	R	–
4.09	Corynebacterium spp.	–	–	–
4.10	Listeria monocytogenes	–	–	–
4.11	Leuconostoc spp., Pediococcus spp.	–	–	–
4.12	Lactobacillus spp. (some species)	–	–	–
4.13	Clostridium ramosum, Clostridium innocuum	–	–	–

Rule no.	Organisms	Quinupristin dalfopristin	Vancomycin	Teicoplanin
4.01	Staphylococcus saprophyticus	–	–	–
4.02	Staphylococcus cohnii,	–	–	–
	Staphylococcus xylosus			
4.03	Staphylococcus capitis	–	–	–
4.04	Other coagulase-negative staphylococci and Staphylococcus aureus	–	–	–
4.05	Streptococcus spp.	–	–	–
4.06	Enterococcus faecalis	R	–	–
4.07	Enterococcus gallinarum, Enterococcus casseliflavus	R	R	–
4.08	Enterococcus faecium	–	–	–
4.09	Corynebacterium spp.	–	–	–
4.10	Listeria monocytogenes	–	–	–
4.11	Leuconostoc spp., Pediococcus spp.	–	R	R
4.12	Lactobacillus spp. (some species)	–	R	R
4.13	Clostridium ramosum, Clostridium innocuum	–	R	–

Rule no.	Organisms	Fosfomycin	Novobiocin	Sulphonamides
4.01	Staphylococcus saprophyticus	R	R	–
4.02	Staphylococcus cohnii,	–	R	–
	Staphylococcus xylosus			
4.03	Staphylococcus capitis	R	–	–
4.04	Other coagulase-negative staphylococci and Staphylococcus aureus	–	–	–
4.05	Streptococcus spp.	–	–	–
4.06	Enterococcus faecalis	–	–	R
4.07	Enterococcus gallinarum, Enterococcus casseliflavus	–	–	R
4.08	Enterococcus faecium	–	–	R
4.09	Corynebacterium spp.	R	–	–
4.10	Listeria monocytogenes	–	–	–
4.11	Leuconostoc spp., Pediococcus spp.	–	–	–
4.12	Lactobacillus spp. (some species)	–	–	–
4.13	Clostridium ramosum, Clostridium innocuum	–	–	–

R. resistant.

a. Low-level resistance to aminoglycosides. Combinations of aminoglycosides with cell wall inhibitors (penicillins and glycopeptides) are synergistic and bactericidal against isolates that are susceptible to cell wall inhibitors and do not display high-level resistance to aminoglycosides.

b. In addition to low-level resistance to aminoglycosides, *Enterococcus faecium* produces a chromosomal AAC(6') enzyme that is responsible for the loss of synergism between aminoglycosides (except gentamicin, amikacin, arbekacin, and streptomycin) and penicillins or glycopeptides.

ACCREDITED FOR ORGANIZATION BY
JOINT COMMISSION INTERNATIONAL

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